CombiScreen® DIP Check

Catalog No. 93010 2 x 15 mL

Lot No.
Expiry

Y 373
2020/05

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IN-VITRO DIAGNOSTIKUM

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Intended Use

FOR IN VITRO DIAGNOSTIC USE

CombiScreen® DIP Check is intended for use as an assayed quality control material for various urinalysis reagent strips and qualitative hCG methods.

GB

Summary and Explanation

The use of quality control materials to objectively monitor the accuracy and precision of procedures is well established in clinical laboratories. CombiScreen® DIP Check is provided at two levels to assist in the monitoring of analytical systems within the clinical range.

Product Description

CombiScreen® DIP Check is human urine-based containing constituents of human and animal origin as well as purified chemicals. Preservatives and stabilizers have been added to maintain product integrity. CombiScreen® DIP Check is a ready-to-use liquid control requiring no reconstitution.

Warnings and Precautions

For in vitro diagnostic use.

Human source material. Treat as potentially infectious.

Each donor unit used in the manufacture of this product has been tested by FDA accepted methods and found non-reactive for the presence of HBsAg and antibody to HIV-1/2, HCV and HIV-1 Ag. While these methods are highly accurate, they do not guarantee that all infected units will be detected. Because no known test method can offer complete assurance the hepatitis B virus, hepatitis C virus, human immunodeficiency virus (HIV) or other infectious agents are absent, all products containing human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.

This product contains 0.09% sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form potentially explosive compounds. Flush with excess water upon disposal.

The material safety data sheet contains further safety-related information. It is available for download from our homepage <http://www.analyticon-diagnostics.com>.

Indications of deterioration

If the controls are turbid or if any evidence of microbial growth or contamination is present, the controls should be discarded. This should be done according to the local guidelines in the same manner as for other biological specimens.

Storage and Stability

CombiScreen® DIP Check is stable until the expiration date on the tube label when stored unopened at 2–8°C. Do not freeze.

Once opened, tubes of control are stable for 75 days when stored tightly capped at 2–8°C or 20 dipstick immersions have taken place, whichever occurs first. CombiScreen® DIP Check may be stored for 30 days at 20–25°C.

Procedure

CombiScreen® DIP Check should be treated the same as patient specimens and run in accordance with the instructions accompanying the test system (instrument, kit, or reagent) being used. Allow the product to reach room temperature prior to use. QC materials should be used in accordance with local, state, and/or federal regulations or accreditation requirements.

Instructions

- Gently invert the tube several times before sampling to ensure homogeneity. This is important to obtain reproducible results.
- For dipstick testing: remove the cap and dip the reagent strip in the tube as if it were a patient sample. In accordance with the manufacturer's instructions, read the reagent strips either visually or with an instrument reader.
For hCG testing: In accordance with the hCG test kit manufacturer's instructions, use controls as if they were patient specimens. For the control of confirmatory pregnancy test kits, be sure to use the transfer pipette included in the kit to deliver the correct sample amount to test device.
- Replace the cap immediately after use and store appropriately.

Materials provided

1x Control L1, 15 ml

1x Control L2, 15 ml

1x Kit insert

Limitations of the Procedure

CombiScreen® DIP Check should not be used past the expiration date on the tube label. CombiScreen® DIP Check is a stabilized liquid product. To obtain consistent assay values, CombiScreen® DIP Check requires storage and handling as detailed in *Storage and Stability*.

Accurate and reproducible results are dependent upon properly functioning instruments and reagents. The ranges given are intended only as a guideline. Each laboratory should establish their own acceptable ranges and tolerance limits based on their test system.

Expected readings

The expected values printed in this insert were derived from replicate analyses of representative samples of the product and are specific to this lot of CombiScreen® DIP Check. Testing data used to establish the expected values were derived from multiple instruments. All values have been assigned with manufacturer's reagents available at the time of assay. Subsequent instrument or reagent modifications may invalidate these expected values. To request a faxed or emailed copy of the value assignment, contact your local distributor or info@analyticon-diagnostics.com.

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Analyticon® Biotechnologies AG