# QuantScopics® Urinalysis Microscopics Control Level 1 (Low) / Level 2 (High)

**LOT** 48301  
48302

**REF** 1481-41  
1482-41

**2019-07**

**IVD**

2°C – 8°C

**LOT**

**IVD**

**Cont.**

**REF**

**EC REP**

European Conformity  
CE-Konformitätsbescheinigung  
Conformité aux normes européennes  
Conformità europea  
Conformidad europea

Lot Number  
Designation du lot  
Numero di lotto  
Denominación de lote

Manufactured by  
Hergestellt von  
Fabriqué par  
Fabricado por

For in vitro diagnostic use  
In-vitro Diagnosticum  
Pour diagnostic in vitro  
Per uso diagnostico in vitro  
Di uso diagnostico in vitro

Biological Risk  
Biosigiführung  
Risque biologique  
Riesgo biológico

Contents of kit  
Inhalt der Packung  
Contenu du coffret  
Contenido della confezione  
Contenido del estuche

Catalog No.  
Bestellf.  
N° de catalogue  
Catalogo n.  
N° de catálogo

Caution, See Product Insert  
Achtung, Siehe Packungsbeilage  
Attention, voir notice d'utilisation  
Attenzione, vedere il foglio  
illustrato del prodotto  
Illustration del prodotto  
Atención, consulte el folleto  
del producto

Authorized Representative  
Bevollmächtigter  
Représentant agréé  
Rappresentante autorizzato  
Representante autorizado

Temperature limitation  
Temperaturbegrenzungen  
Limites de température  
Limiti di temperatura  
límite de temperatura

Consult instructions for use  
Gebrauchsanweisung beachten  
Consultez les instructions d'utilisation  
Consultare le istruzioni d'uso  
Consulte las instrucciones de uso

Use by (last day of month)  
Verwendbar bis letzter Tag des Monats  
Utilisable jusqu'à (dernier jour du mois indiqué)  
Da utilizzare prima del (ultimo giorno del mese)  
Estable hasta (ultimo día del mes)

## English

### Intended use

The QuantScopics Urinalysis Microscopics Control is intended for use as quality control material for the microscopic evaluation of various urine sediment methods.

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method. Microscopic QC controls should be run each day the test is performed.<sup>1</sup>

### Summary and Explanation

Standardized microscopic evaluation of urine sediment is an important part of routine analysis of urine. Along with physical and chemical analysis, microscopic examination of urine can provide valuable information regarding not only renal and urinary tract disease, but also metabolic diseases unrelated to the kidney.<sup>2</sup> Urinary sediment microscopy generally includes the detection and identification of red blood cells, leukocytes, epithelial cells, bacteria, casts and crystals.<sup>3,4</sup>

### Product Description

The QuantScopics Urinalysis Microscopics Control is supplied in two levels, 4 x 120 mL bottles per level per box. It is liquid, ready-to-use, requiring no reconstitution or dilution. The control is made from human urine to which stabilized human red blood cells, human white blood cells, and uric acid, calcium oxalate and calcium phosphate crystals have been added. Preservatives have been added to maintain product integrity.

### Warnings and Precautions

**POTENTIAL BIOHAZARDOUS MATERIAL.** All blood donor units comprising the human source material used in the manufacture of this product have been tested and found non-reactive for Hepatitis B surface antigen, Hepatitis C, and HIV antibody when tested by FDA accepted methods. No known test method can assure that a product derived from human blood does not contain Hepatitis or HIV virus. It is recommended that such samples be handled according to the Centers for Disease Control's Bio-Safety Level 2 recommendations.

### Warning ⚠ Hazard (H) and Precautionary (P) Statements

Contains Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone, Benzyl alcohol, Chloramphenicol, Novobiocin, Erythromycin.

H317 – May cause an allergic skin reaction

P261 – Avoid breathing vapors, mist, or spray.

P280 – Wear protective gloves, protective clothing, and eye protection.

P302+P352 – IF ON SKIN: Wash with plenty of water.

P321 – Specific Treatment (see section 4 on the SDS).

P333+P313 – If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 – Take off contaminated clothing and wash it before reuse.

P501 – Dispose of contents/container in accordance with local, regional, national, and international regulations.

Safety Data Sheet (SDS) available for professional users at [quantimetrix.com](http://quantimetrix.com)

### Storage and Stability

The QuantScopics Urinalysis Microscopics Control should be stored tightly capped at 2°C – 8°C when not in use. Do not freeze. When stored at 2°C – 8°C the controls are stable until the expiration date stated on the label. Discard the controls if turbid or any evidence of microbial contamination is present.

### Procedure

Remove the control from the refrigerator and replace the cap on the control bottle with the pour seal cap included in the control box. Allow the control to come to room temperature (18°C – 25°C) for approximately 15 – 90 minutes depending on the volume remaining in the bottle. Mix the control thoroughly by inverting the bottle at least 20 times to assure homogeneity of the contents. Avoid foaming. Thorough mixing with each use is important in order to obtain reproducible results. Treat the control as you would a patient sample in accordance with the manufacturer's instructions for the standardized microscopic urinalysis system you are using. The National Committee for Clinical Laboratory Standards (NCCLS) recommends the use of standardized systems in order to yield standardized, reproducible results and to enable the reporting of abnormal sediment elements per unit volume.<sup>1</sup> Immediately close the pour seal cap and store the control at 2°C – 8°C when not in use.

### Expected Ranges

The expected ranges for each type of formed element found in this lot of QuantScopics Control is listed in the table on page 4. These ranges were determined by assay of multiple bottles of the indicated lot by the methods listed. A 12 mL sample volume was used for each testing method. The samples were centrifuged at 400 RCF (relative centrifugal force) for 5 minutes. After centrifugation, urine sediment was resuspended in either ~0.5 or ~1.0 mL of remaining supernatant according to the respective cell count method procedure. The ranges listed are based on the range of elements observed in 10 high power fields. Use of other systems or protocols may yield differing results. Each laboratory should establish its own precision parameters.

### Limitations

Any future changes made by the manufacturer of a test system may give different values from the indicated range. Detailed information on the limitations of each test system is included in the limitations section of the system manufacturer's package insert. Technical updates can be found on our

website. The Quality Control Log will no longer be shipped with the product. You can obtain it from the Quantimetrix website at [www.quantimetrix.com](http://www.quantimetrix.com) or contact Tech Support at (310) 536-0006, option 3.

### References

1) NCCLS Document GP16-A. Routine urinalysis and collection, transportation, and preservation of urine specimens; approved guideline. Wayne PA: National Committee for Clinical Laboratory Standards, 1995. 2) Bradley GM, Benson ES. Todd-Sanford Clinical Diagnosis by Laboratory Methods, 15th ed. Philadelphia PA: Saunders, 1974. 3) Bologna CV. Understanding Laboratory Medicine. St. Louis MO: CV Mosby Co, 1971. 4) Dudas HC. Lab Med 1981;12:765.

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Gemisch, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone, Benzyl alcohol, Chloramphenicol, Novobiocin, Erythromycin.

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# QuantScopics® Urinalysis Microscopics Control

## Level 1 (Low) / Level 2 (High)

Microscopic Evaluation of Urine Sediment	Mikroskopische Beurteilung von Urinsediment	Evaluation microscopique du sédiment urinaire	Valutazione microscopica del sedimento urinario	Evaluación microscópica de la sedimentación presente en la orina
Analytes / Method		Level 1 48301	Level 2 48302	
<b>Red Blood Cells</b>				
UriSed Variants / COBIO Variants		0 - 30/µL	20 - 240/µL	
IRIS® Diagnostics iQ200 Series Analyzers		1 - 16/µL	35 - 98/µL	
AVE-764B Automated Urine Formed Elements Analyzer		0 - 19/µL	45 - 118/µL	
Kova® Glassic™ Slide 10 with grids		1 - 13/µL	16 - 89/µL	
Non-grid slides (~0.5 mL resuspended sediment)		0 - 6/hpf	1 - 29/hpf	
Non-grid slides (~1.0 mL resuspended sediment)		0 - 6/hpf	2 - 29/hpf	
Cen-Slide® System		0 - 12/hpf	2 - 35/hpf	
Slide & Coverslip (~0.5 mL resuspended sediment)		0 - 6/hpf	2 - 22/hpf	
Slide & Coverslip (~1.0 mL resuspended sediment)		0 - 6/hpf	1 - 20/hpf	
FisherBrand UriSystem DeciSlide		0 - 6/hpf	0 - 20/hpf	
<b>White Blood Cells</b>				
UriSed Variants / COBIO Variants		0 - 20/µL	20 - 100/µL	
IRIS Diagnostics iQ200 Series Analyzers		2 - 14/µL	38 - 91/µL	
AVE-764B Automated Urine Formed Elements Analyzer		0 - 10/µL	20 - 70/µL	
Kova Glassic Slide 10 with grids		0 - 9/µL	6 - 75/µL	
Non-grid slides (~0.5 mL resuspended sediment)		0 - 10/hpf	0 - 31/hpf	
Non-grid slides (~1.0 mL resuspended sediment)		0 - 5/hpf	0 - 19/hpf	
Cen-Slide System		0 - 11/hpf	2 - 25/hpf	
Slide & Coverslip (~0.5 mL resuspended sediment)		0 - 5/hpf	0 - 16/hpf	
Slide & Coverslip (~1.0 mL resuspended sediment)		0 - 6/hpf	0 - 14/hpf	
FisherBrand UriSystem DeciSlide		0 - 4/hpf	0 - 18/hpf	
<b>Casts</b>				
UriSed Variants / COBIO Variants		none	none	
IRIS Diagnostics iQ200 Series Analyzers		none	none	
AVE-764B Automated Urine Formed Elements Analyzer		none	none	
Kova Glassic Slide 10 with grids		none	none	
Non-grid slides (~0.5 mL resuspended sediment)		none	none	
Non-grid slides (~1.0 mL resuspended sediment)		none	none	
Cen-Slide System		none	none	
Slide & Coverslip (~0.5 mL resuspended sediment)		none	none	
Slide & Coverslip (~1.0 mL resuspended sediment)		none	none	
<b>Crystals</b>				
UriSed Variants / COBIO Variants		may be present	may be present	
IRIS Diagnostics iQ200 Series Analyzers		Present	Present	
AVE-764B Automated Urine Formed Elements Analyzer		may be present	may be present	
Kova Glassic Slide 10 with grids		Present	Present	
Non-grid slides (~0.5 mL resuspended sediment)		Present	Present	
Non-grid slides (~1.0 mL resuspended sediment)		Present	Present	
Cen-Slide System		Present	Present	
Slide & Coverslip (~0.5 mL resuspended sediment)		Present	Present	
Slide & Coverslip (~1.0 mL resuspended sediment)		Present	Present	
<b>Bacteria</b>				
UriSed Variants / COBIO Variants		may be present	may be present	
IRIS Diagnostics iQ200 Series Analyzers		may be present	may be present	
AVE-764B Automated Urine Formed Elements Analyzer		may be present	may be present	
Kova Glassic Slide 10 with grids		may be present	may be present	
Non-grid slides (~0.5 mL resuspended sediment)		may be present	may be present	
Non-grid slides (~1.0 mL resuspended sediment)		may be present	may be present	
Cen-Slide System		may be present	may be present	
Slide & Coverslip (~0.5 mL resuspended sediment)		may be present	may be present	
Slide & Coverslip (~1.0 mL resuspended sediment)		may be present	may be present	

Quantimetrix Corporation  
2005 Manhattan Beach Blvd.  
Redondo Beach, CA 90278  
+1.310.536.0006

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MDSS GmbH