Ammonia Alcohol Control / Level 1, Level 2, Level 3

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English

Intended Use

The Quantimetrix Ammonia/Alcohol control is intended as a means of monitoring various ammonia and ethanol assay methods to validate quantitation of patient samples.

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method.

Product Description

The Quantimetrix Ammonia/Alcohol Control is supplied in three levels, 3 x 5 mL, as a ready-to-use liquid requiring no reconstitution or dilution. It is fortified with ethanol and reagent grade chemicals. Preservatives have been added to inhibit microbial growth.

Precautions

For in vitro diagnostic use only. Observe the same safety precautions used when handling potentially infectious biological material.

Storage and Stability

Store the control at 2°C–8°C. When stored unopened at 2°C–8°C, the control is stable until the expiration date stated on the label. The control has an open vial stability of four months after opening and use. Discard the control if it is turbid or there is any evidence of microbial contamination. Discard controls in the same manner as other biological specimens, according to local guidelines.

Procedure

Remove the controls from the refrigerator and allow to come to room temperature (18°C–25°C), at least 15 minutes, depending on remaining volume. Gently invert the control to assure homogeneity of the contents. Avoid foaming. Treat the control as you would a patient sample in accordance with the manufacturer's requirements of the test method. Immediately recap the control and return it to 2°C–8°C when not in use.

Expected Range

The expected ranges have been established from interlaboratory data. Each laboratory should establish its own precision parameters. A quality control program, Quantrol,[®] is available from Quantimetrix. Please inquire.

Limitations

Any future changes made by the manufacturer of a test method may give different values from the indicated range. Limitations of the test method are included in the package insert for the reagent kit or instrument being used. Technical updates can be found on our website. Safety Data Sheet (SDS) available for professional users at quantimetrix.com.

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2 Siemens, Inc., Dimension® ExL
Quantimetrix Corporation
MDSS GmbH