
Quantimetrix® DipperPOCT® Urinalysis Dipstick Control

1 & 2

Level

Cont. 10 Pouches - Level 1
10 Pouches - Level 2

REF 1400-02 2021-02-28

Quantimetrix
Redondo Beach, CA 90278

EC REP MDSS GmbH, Schiffgraben 41
30175 Hannover, Germany
U.S. Pat. No. 9,835,562

	European Conformity CE-Kennzeichnungsverfahren Conformità ai norme europee Conformité europea Conformidad europea		Caution, See Product Insert Achtung, Siehe Packungsbeilage Attention, voir notice d'utilisation Attenzione, vedere il foglio Illustrativo del prodotto Atención, consulte el folleto del producto
REF	Catalog No. Besteller N° de catalogue Catalogo n. N° de catálogo	IVD	For in vitro diagnostic use In-vitro Diagnostikum Pour diagnostic in vitro Per uso diagnostico in vitro De uso diagnóstico in vitro
	Manufactured by Hergestellt von Fabriqué par Fabricado da Fabricado por		Temperature limitation Temperaturbegrenzungen Limites de température límites de temperatura
EC REP	Authorized Representative Bevollmächtigter Représentant agréé Representante autorizab Representante autorizado	Cont.	Contents of kit Inhalt der Packung Contenu du coffret Contenuto della confezione Contenido del estuche
	Consult instructions for use Gebrauchsanweisung beachten Consulter les instructions d'utilisation Consultare le istruzioni l'uso Consulte las instrucciones de uso		Use by (last day of month) Verwendbar bis (letzter Tag des Monats) Utilisable jusqu'à (dernier jour du mois indiqué) Da utilizzare prima del (ultimo giorno del mese) Estable hasta (último día del mes)
	Single-use only (do not reuse) Nur zum einmaligen Gebrauch (nicht wiederverwenden) À usage unique (ne pas réutiliser) Esclusivamente monouso (non riutilizzare) Un solo uso (no reutilizar)	LOT	Lot Number Bezeichnung Designation du lot Número di lotto Denominación de lote
Rx only	U.S. Federal law restricts this device to sale by or only on the order of a licensed healthcare practitioner. Gemäß US-Bundesgesetz darf dieses Gerät nur von einer lizenzierten Fachkraft im Gesundheitswesen bzw. in deren Auftrag verkauft. La loi fédérale des États-Unis restreint la vente de cet appareil par ou uniquement sur ordonnance d'un professionnel de santé agréé. Legge federale degli Stati Uniti limita la vendita del dispositivo da parte di o solo su prescrizione di medici specialisti. La ley federal de Estados Unidos permite la venta de este dispositivo exclusivamente a médicos o bajo prescripción facultativa.		

English Intended Use

The Quantimetrix Dipper POCT® Urinalysis Dipstick Control is an assayed quality control material used to monitor the performance of urinalysis test strips containing various combinations of reagent pads for bilirubin, blood, creatinine, glucose, ketones, leukocytes, microalbumin, nitrite, pH, protein, specific gravity, and urobilinogen. Dipper POCT may also serve as a control for confirmatory tests such as **K-CHECK**, and **Ictotest®** reagent tablets as well as for **hCG** lateral flow methods.

Summary and Explanation

Control materials having known component concentrations are an integral part of diagnostic procedures. Daily monitoring of control values establishes intra-laboratory parameters for accuracy and precision of the test method.

Product Description

The Dipper POCT is a two level urinalysis dipstick control. Each level is contained in a single-use plastic pouch filled with 1.5mL of liquid control material that is ready-to-use, requiring no reconstitution or dilution.

The controls are formulated from a simulated urine matrix fortified with preservatives and stabilizers to maintain product integrity and inhibit microbial growth. Each control level is formulated to target levels with compounds that produce the desired reaction when tested by the methods indicated in the **Intended Use section**.

Caution

Contains human chorionic gonadotropin (hCG) from pregnancy urine. The human hCG source material used in the manufacture of this product has been tested and found non-reactive for Hepatitis B Surface Antigen, Hepatitis C, and HIV 1&2 antibody when tested with FDA accepted methods. No known test method can assure that a product derived from human material does not contain Hepatitis or HIV virus. Handle the QC material as you would a patient sample. QC materials should be used and disposed of in accordance with regulatory and accreditation requirements.

Warning ⚠ Hazard (H) and Precautionary (P) Statements

Contains Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone.
1,2-Propylene Glycol, Level 1; 2,4-Pentanedione, Level 2.

H317 – May cause an allergic skin reaction

P261 – Avoid breathing vapors, mist, or spray.

P272 – Contaminated work clothing should not be allowed out of the workplace.

P280 – Wear protective gloves, protective clothing, and eye protection.

P302+P352 – IF ON SKIN: Wash with plenty of water.

P333+P313 – If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 – Take off contaminated clothing and wash it before reuse.

P501 – Dispose of contents/container in accordance with local, regional, national, and international regulations.

Safety Data Sheet (SDS) available for professional users at quantimetrix.com

Storage and Stability

The Dipper POCT Urinalysis Dipstick Control kit should be stored at 2°C-8°C and is stable until the expiration date stated on the label. **Do not freeze**. Alternatively, each kit may be stored at room temperature (18°C-25°C) for up to 3 months. Kits and/or loose pouches stored at room temperature should be marked with the 3 month expiration date per each laboratory's standard practices. Once opened, pouches should be used within 30 minutes.

Procedure for Dipstick Testing

- 1 Remove a Dipper POCT pouch from the kit. If the kits have been stored at 2°C-8°C, allow the pouch to equilibrate to room temperature for at least 10 minutes prior to use.
- 2 Remove a urine dipstick from the container and set aside in an easily accessed location.
- 3 Hold the Dipper POCT pouch upright and allow the fluid to drain into the lower portion of the pouch. A gentle tap or shake may be required to move the fluid. Become familiarized with the peelable chevron seal and the STOP-PEEL POINT indicated by the Top of the Ribbed Area. (Figure 1.)
- 4 Use two hands to separate the clear flap from the opaque flap. (Figure 2.)
- 5 While opening, take care not to apply pressure to the fluid channel as the contents will squirt out. (Figure 3.)
- 6 Grab each flap placing your thumb on the inside center of the flaps. Apply gentle and consistent pressure to peel open the pouch until you reach the STOP-PEEL POINT. Once you've reached the STOP-PEEL POINT, you will not be able to peel the package further. A triangular shaped opening along the clear flap should be exposed. (Figure 4.)
- 7 Insert the dipstick and visually verify that all of the reagent pads have been immersed into the fluid. Remove the dipstick, discard the pouch, and follow the manufacturer's instructions for reading the dipstick. (Figure 5.)

Procedure for hCG and Confirmatory Tablet Testing

To use the Dipper POCT control for hCG lateral flow and confirmatory tablet methods, begin by following steps 1-6 above. Once the pouch has been opened, a pipette or similar device may be inserted into the pouch to remove the required volume per the hCG test kit or confirmatory tablet test manufacturer's instructions. Discard each pouch after each usage.

Expected Values

All expected value ranges and assignments are determined using pooled data from the system manufacturer and/or internal data generated as per the manufacturer's instructions.

Limitations

Each pouch of Dipper POCT is single use only and should be discarded after each use. Once the control fluid is removed for hCG or confirmatory tablet testing, that pouch should not be used for any further testing. Similarly, once the pouch is used for a urine test strip, it should not be used for any further testing.

Any future changes made by the manufacturer of a test method may give different values from the indicated range. Detailed information on the limitations of each test method is included in the limitations section of the manufacturers' package insert. Technical updates can be found on our website. The Quality Control Log will no longer be shipped with the product. You can obtain it from the Quantimetrix website at quantimetrix.com or contact Tech Support at (310) 536-0006, option 3.

Siemens Multistix Users

Colors produced by the **urobilinogen** and/or **bilirubin** reactions on these dipsticks with the urinalysis dipstick control may not be characteristic of those shown on the manufacturer's label when reading the dipstick reactions visually. The urobilinogen reactions are consistent and intensity with the increase in the urobilinogen concentration but may not provide an exact color match to those displayed on the label.

Note: Siemens® CLINITEK 50 and Siemens® STATUS or CLINITEK STATUS PLUS may see an Albumin/Creatinine ratio result of "Abnormal" with the Level 1 control.

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Gemisch, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone.

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IMPORTANT: FOLLOW THIS PROCEDURE FOR DIPSTICK TESTING

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
Figure 4 / Abbildung 4 / Figura 4


Figure 5 / Abbildung 5 / Figura 5

Analytes	Level 1 - Lot 40031	Level 2 - Lot 40032
Confirmatory and Other Tests		
K-CHECK (Ketones)	Negative	Small - Large
Ictotest (Bilirubin)	Negative	Positive
Refractometer (Specific Gravity)	1.020 - 1.027	1.013 - 1.020
hCG	Negative	Positive
pH Paper	4 - 6	7 - 9
Sulfosalicylic Acid (Total Protein) ⁵	Negative (≤ 0.0511)	Positive (≥ 0.5011)
Siemens® VISUAL TESTING (Visual Test Strips Only)		
Glucose	Negative	100 - 1000 mg/dL (Tr - 3+)
Bilirubin	Negative	Small - Large (1+ - 3+)
Ketones	Negative	15 - 160 mg/dL (Sm - Lg)
Specific Gravity	1.020 - 1.030	1.005 - 1.020
Blood	Negative	Small - Large (1+ - 3+)
pH	5 - 6.5	7 - 9
Protein	Negative	30 - 300 mg/dL (1+ - 3+)
Urobilinogen	Normal (0.2 E.U./dL) ⁴	1.0 - 8.0 E.U./dL*
Nitrite	Negative	Positive
Leukocytes	Negative	Trace - Large (Tr - 3+)
Creatinine ²	10 - 100 mg/dL	100 - 300 mg/dL
Siemens® CLINITEK 50		
Glucose	Negative	100 - ≥ 1000 mg/dL (Tr - 3+)
Bilirubin	Negative	Small - Large (1+ - 3+)
Ketones	Negative	15 - ≥ 80 mg/dL (1+ - 3+)
Specific Gravity	1.020 - ≥ 1.030	≤ 1.005 - 1.020
Blood	Negative	Small - Large (1+ - 3+)
pH	5 - 6.5	7 - 8.5
Protein	Negative	30 - ≥ 300 mg/dL (1+ - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	1.0 - 4.0 E.U./dL
Nitrite	Negative	Positive
Leukocytes	Negative	Trace - Large (Tr - 3+)
Microalbumin ¹	≤ 10 - 30 mg/L	30 - 300 mg/L
Creatinine ³	10 - 100 mg/dL	100 - 300 mg/dL
Siemens® CLINITEK 500		
Glucose	Negative	100 - ≥ 1000 mg/dL (Tr - 3+)
Bilirubin	Negative	Small - Large (1+ - 3+)
Ketones	Negative	15 - ≥ 80 mg/dL (1+ - 3+)
Specific Gravity	1.020 - ≥ 1.030	≤ 1.005 - 1.020
Blood	Negative	Small - Large (1+ - 3+)
pH	5 - 6.5	7 - ≥ 9
Protein	Negative	30 - ≥ 300 mg/dL (1+ - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	1.0 - ≥ 8.0 E.U./dL
Nitrite	Negative	Positive
Leukocytes	Negative	Trace to Large (Tr - 3+)
Creatinine ²	10 - 100 mg/dL	100 - 300 mg/dL
Siemens® CLINITEK ADVANTUS		
Glucose	Negative	100 - ≥ 1000 mg/dL (Tr - 3+)
Bilirubin	Negative	Small - Large (1+ - 3+)
Ketones	Negative	15 - ≥ 80 mg/dL (1+ - 3+)
Specific Gravity	1.020 - ≥ 1.030	≤ 1.005 - 1.020
Blood	Negative	Small - Large (1+ - 3+)
pH	5 - 6.5	7 - ≥ 9
Protein	Negative	30 - ≥ 300 mg/dL (1+ - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	1.0 - ≥ 8.0 E.U./dL
Nitrite	Negative	Positive
Leukocytes	Negative	Trace to Large (Tr - 3+)
Creatinine ²	10 - 100 mg/dL	100 - 300 mg/dL
Siemens® CLINITEK STATUS or Siemens® CLINITEK STATUS PLUS		
Glucose	Negative	100 - ≥ 1000 mg/dL (Tr - 3+)
Bilirubin	Negative	Small to Large (1+ - 3+)
Ketones	Negative	15 - ≥ 160 mg/dL (1+ - 4+)
Specific Gravity	1.020 - ≥ 1.030	1.005 - 1.025
Blood	Negative	Small to Large (1+ - 3+)
pH	5 - 6.5	7 - ≥ 9
Protein	Negative	30 - ≥ 300 mg/dL (1+ - 3+)
Urobilinogen	Normal (0.2 E.U./dL)	1.0 - ≥ 8.0 E.U./dL
Nitrite	Negative	Positive
Leukocytes	Negative	Trace to Large (Tr - 3+)
Microalbumin ¹	10 - 30 mg/L	30 - 300 mg/L
Creatinine ³	10 - 100 mg/dL	100 - 300 mg/dL
hCG	Negative	Positive

Footnotes for values

- 1 Values only apply to Clinitek Microalbumin Reagent Strips when read on the Clinitek 50 and Status.
 - 2 Values only apply to Multistix Pro Reagent Strips.
 - 3 Values only apply to Multistix Pro and Clinitek Microalbumin Reagent Strips when read on Clinitek Urine Analyzers.
 - 4 Some customers may obtain false positives.
 - 5 Absorbance at 620 nm.
- * See limitations.

Level 1

LOT 40031

Level 2

LOT 40032

Level 1&2 Expiration Date

2021-02-28

Representative of dispensing method.
Vertreter der Verpackungsart.
Représentant de la méthode de distribution.
Rappresentante del metodo di erogazione.
Representante del método de dispensación.

Dipper and Dipper POCT
are registered trademarks
of Quantimetrix.

Quantimetrix Corporation
2005 Manhattan Beach Blvd.
Redondo Beach, CA 90278
+1.310.536.0006

IE-MO39030A-05/18