Telefax: +49 (0) 6454/7991-30

## **Safety Data Sheet**

according to Regulation (EC) No 1907/2006

# CombiScreen® Dip Check Level 1

Revision date: 17.02.2023 Product code: 1R93010 Page 1 of 7

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

CombiScreen® Dip Check Level 1

Product group: Endprodukt / Endproduct

## 1.2. Relevant identified uses of the substance or mixture and uses advised against

### 1.3. Details of the supplier of the safety data sheet

Company name: Analyticon® Biotechnologies GmbH

Street: Am Mühlenberg 10
Place: D-35104 Lichtenfels
Telephone: +49 (0) 6454/7991-0

E-mail: support@analyticon-diagnostics.com

Contact person: Zentrale Telephone: +49 (0) 6454/7991-0

Internet: http://www.analyticon-diagnostics.com

**1.4. Emergency telephone** Zentrale: +49 (0) 6454/7991-0

number:

#### **SECTION 2: Hazards identification**

## 2.1. Classification of the substance or mixture

### Regulation (EC) No 1272/2008

This mixture is not classified as hazardous in accordance with Regulation (EC) No 1272/2008.

### 2.2. Label elements

## Regulation (EC) No 1272/2008

# Special labelling of certain mixtures

Restricted to professional users.

### **SECTION 3: Composition/information on ingredients**

### 3.2. Mixtures

## **Hazardous components**

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification (Regulation (EC) No 1272/2008)			
	Human Source Material			10-60 %

Full text of H and EUH statements: see section 16.

# **Further Information**

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].

## **SECTION 4: First aid measures**

## 4.1. Description of first aid measures

### After inhalation

Provide fresh air.

### After contact with skin

Wash with plenty of water. Take off contaminated clothing and wash it before reuse.

according to Regulation (EC) No 1907/2006

## CombiScreen® Dip Check Level 1

Revision date: 17.02.2023 Product code: 1R93010 Page 2 of 7

### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

#### After ingestion

Rinse mouth immediately and drink plenty of water.

Rinse mouth thoroughly with water.

Seek medical advice immediately.

# 4.2. Most important symptoms and effects, both acute and delayed

No information available.

#### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

### **SECTION 5: Firefighting measures**

### 5.1. Extinguishing media

### Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

The product itself does not burn.

## 5.2. Special hazards arising from the substance or mixture

Non-flammable.

#### 5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

### **SECTION 6: Accidental release measures**

## 6.1. Personal precautions, protective equipment and emergency procedures

### General advice

Wear breathing apparatus if exposed to vapours/dusts/aerosols.

#### 6.2. Environmental precautions

Do not allow to enter into surface water or drains.

Prevent spread over a wide area (e.g. by containment or oil barriers).

Do not allow to enter into soil/subsoil.

## 6.3. Methods and material for containment and cleaning up

### Other information

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

Take up dust-free and set down dust-free.

### 6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## **SECTION 7: Handling and storage**

## 7.1. Precautions for safe handling

### Advice on safe handling

Use only in well-ventilated areas.

The floor should be leak tight, jointless and not absorbent.

All work processes must always be designed so that the following is excluded:

## Advice on protection against fire and explosion

Usual measures for fire prevention.

When using do not smoke.

### Advice on general occupational hygiene

Take off contaminated clothing. Wash hands before breaks and after work. When using do not eat or drink.

according to Regulation (EC) No 1907/2006

# CombiScreen® Dip Check Level 1

Revision date: 17.02.2023 Product code: 1R93010 Page 3 of 7

Wash hands before breaks and after work.

### Further information on handling

When using do not eat, drink, smoke, sniff.

## 7.2. Conditions for safe storage, including any incompatibilities

### Requirements for storage rooms and vessels

Keep container tightly closed.

Keep only in the original container in a cool, well-ventilated place.

#### Further information on storage conditions

2 8

Protect against:

# **SECTION 8: Exposure controls/personal protection**

### 8.1. Control parameters

### 8.2. Exposure controls



### Individual protection measures, such as personal protective equipment

### Eye/face protection

Wear eye/face protection.

#### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. EN ISO 374

Breakthrough times and swelling properties of the material must be taken into consideration.

### Skin protection

Wear suitable protective clothing. Suitable protective clothing:

#### Respiratory protection

In case of inadequate ventilation wear respiratory protection. Respiratory protection necessary at:

# **SECTION 9: Physical and chemical properties**

# 9.1. Information on basic physical and chemical properties

Physical state:

Colour:

Odour: characteristic

Test method

Melting point/freezing point:

Boiling point or initial boiling point and

not determined

100 °C

boiling range:

Print date: 31.05.2023

## **Safety Data Sheet**

according to Regulation (EC) No 1907/2006

## CombiScreen® Dip Check Level 1

Revision date: 17.02.2023 Product code: 1R93010 Page 4 of 7

Flammability:

not determined
not applicable

Decomposition temperature:

not determined

Partition coefficient n-octanol/water:

Vapour pressure:

not determined

Density:

not determined

Relative vapour density:

not determined

#### 9.2. Other information

## Information with regard to physical hazard classes

Explosive properties

The study does not need to be conducted because there are no chemical groups associated with explosive properties present in the molecule.

Self-ignition temperature

Solid: not determined
Gas: not applicable

Other safety characteristics

Evaporation rate: not determined Solid content: not determined

## **SECTION 10: Stability and reactivity**

## 10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

#### 10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

## 10.3. Possibility of hazardous reactions

No known hazardous reactions.

## 10.4. Conditions to avoid

May cause decomposition by long-term light influence.

#### 10.5. Incompatible materials

No information available.

## **SECTION 11: Toxicological information**

## 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

#### Acute toxicity

No information available.

## **ATEmix calculated**

ATE (oral) > 2000 mg/kg; ATE (dermal) > 2000 mg/kg; ATE (inhalation vapour) > 20 mg/l; ATE (inhalation dust/mist) > 5 mg/l

### Sensitising effects

No information available.

### STOT-repeated exposure

No information available.

## Specific effects in experiment on an animal

No information available.

#### Additional information on tests

The mixture is classified as not hazardous according to Directive 1999/45/EC.

## 11.2. Information on other hazards

according to Regulation (EC) No 1907/2006

## CombiScreen® Dip Check Level 1

Revision date: 17.02.2023 Product code: 1R93010 Page 5 of 7

### **Endocrine disrupting properties**

No information available.

# **SECTION 12: Ecological information**

## 12.1. Toxicity

The product has not been tested.

### 12.2. Persistence and degradability

No data available

### 12.3. Bioaccumulative potential

No data available

### 12.4. Mobility in soil

No data available

#### 12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.

The product has not been tested.

#### 12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

### 12.7. Other adverse effects

No data available

#### **Further information**

Avoid release to the environment.

### **SECTION 13: Disposal considerations**

### 13.1. Waste treatment methods

#### **Disposal recommendations**

Dispose of waste according to applicable legislation. Dispose of waste according to applicable legislation.

Dispose of waste according to "Kreislaufwirtschafts- und Abfallgesetz (KrW-/AbfG)".

# List of Wastes Code - residues/unused products

160506 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and

discarded chemicals; laboratory chemicals, consisting of or containing hazardous substances,

including mixtures of laboratory chemicals; hazardous waste

## List of Wastes Code - used product

160506 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and

discarded chemicals; laboratory chemicals, consisting of or containing hazardous substances,

including mixtures of laboratory chemicals; hazardous waste

# List of Wastes Code - contaminated packaging

150106 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND

PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately

collected municipal packaging waste); mixed packaging

### Contaminated packaging

Wash with plenty of water. Completely emptied packages can be recycled.

# **SECTION 14: Transport information**

### Land transport (ADR/RID)

14.1. UN number or ID number:No dangerous good in sense of this transport regulation.14.2. UN proper shipping name:No dangerous good in sense of this transport regulation.14.3. Transport hazard class(es):No dangerous good in sense of this transport regulation.14.4. Packing group:No dangerous good in sense of this transport regulation.

according to Regulation (EC) No 1907/2006

## CombiScreen® Dip Check Level 1

Revision date: 17.02.2023 Product code: 1R93010 Page 6 of 7

Inland waterways transport (ADN)

14.1. UN number or ID number:No dangerous good in sense of this transport regulation.14.2. UN proper shipping name:No dangerous good in sense of this transport regulation.14.3. Transport hazard class(es):No dangerous good in sense of this transport regulation.14.4. Packing group:No dangerous good in sense of this transport regulation.

Marine transport (IMDG)

14.1. UN number or ID number:No dangerous good in sense of this transport regulation.14.2. UN proper shipping name:No dangerous good in sense of this transport regulation.14.3. Transport hazard class(es):No dangerous good in sense of this transport regulation.14.4. Packing group:No dangerous good in sense of this transport regulation.

Air transport (ICAO-TI/IATA-DGR)

14.1. UN number or ID number:No dangerous good in sense of this transport regulation.14.2. UN proper shipping name:No dangerous good in sense of this transport regulation.14.3. Transport hazard class(es):No dangerous good in sense of this transport regulation.14.4. Packing group:No dangerous good in sense of this transport regulation.

14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

14.6. Special precautions for user

No dangerous good in sense of this transport regulation.

#### 14.7. Maritime transport in bulk according to IMO instruments

No dangerous good in sense of this transport regulation.

## **SECTION 15: Regulatory information**

## 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

## National regulatory information

Water hazard class (D): -- non-hazardous to water

## 15.2. Chemical safety assessment

Chemical safety assessments for substances in this mixture were not carried out.

## **SECTION 16: Other information**

### Changes

This data sheet contains changes from the previous version in section(s): 11.

### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route

(European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

## **Further Information**

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations. The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing,

according to Regulation (EC) No 1907/2006

# CombiScreen® Dip Check Level 1

Revision date: 17.02.2023 Product code: 1R93010 Page 7 of 7

transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights.

The receiver of our product is singularly responsible for adhering to existing laws and regulations.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)