

CombiScreen® Dip Check Controls

Urinalysis Dipstick and hCG Control



Ready to Use Dipper Control

CombiScreen® Dip Check is a reliable quality control based on human standard material, for both visual and instrument read urine dipsticks. The control is very easy-to-use as it simulates patient testing with the full-immersion method, and come in shatter-proof vials with 2D barcodes for direct entry of lot number.

Includes values for all dipstick analytes plus Microalbumin and Creatinine, and qualitative results for hCG early pregnancy detection test methods.

- Open vial stability: 20 dips
- 75 days stability at 2–8° C
- 30 days stability at room temperature

Analytes

Bilirubin	Ketones	Specific Gravity
Blood	Leukocytes	Urobilinogen
Creatinine	Nitrite	Nitrite
Glucose	pH	Creatinine
hCG	Protein	



Learn More

[Combiscreen® Dip Check Webpage](#)

Contact

Email sales@sterilab.co.uk
Telephone 01423 523300
Website www.sterilab.co.uk

