



Type	Reflectance photometer with 4 wavelengths: 505, 530, 620, 660 nm
Throughput	Up to 500 tests/hour
Memory	Patient test results: 5000 / QC test results: 1000
Display	5.7" LCD VGA capacitive touchscreen (resolution: 640x480)
Interfaces	PS2 (external keyboard, barcode reader), serial RS232, USB Type B, USB Type A, Ethernet, optional external Wi-Fi connection
Printer	Internal thermal printer (roll diameter max. 60 mm)
Communication protocols	LIS2 (ASTM+), HL7
Dimensions	300 x 300 x 170 mm (WxDxH)
Weight	6 kg
Power supply	IN: 100–240 V AC, 50/60 Hz OUT: 12 V DC, 5 A